

---

# **franka\_ros\_interface**

## **Documentation**

**Release 0.0.1**

**Saif Sidhik**

**Mar 16, 2020**



## Contents:

<b>1</b>	<b>Python API Documentation</b>	<b>1</b>
1.1	franka_interface . . . . .	1
1.1.1	ArmInterface . . . . .	1
1.1.2	GripperInterface . . . . .	6
1.1.3	RobotEnable . . . . .	9
1.1.4	RobotParams . . . . .	9
1.2	franka_moveit . . . . .	10
1.2.1	PandaMoveGroupInterface . . . . .	10
1.2.2	ExtendedPlanningSceneInterface . . . . .	12
1.3	franka_tools . . . . .	12
1.3.1	CollisionBehaviourInterface . . . . .	12
1.3.2	FrankaControllerManagerInterface . . . . .	13
1.3.3	ControllerParamConfigClient . . . . .	16
1.3.4	FrankaFramesInterface . . . . .	18
1.3.5	JointTrajectoryActionClient . . . . .	19
<b>2</b>	<b>Indices and tables</b>	<b>21</b>
	<b>Python Module Index</b>	<b>23</b>
	<b>Index</b>	<b>25</b>



## Python API Documentation

### 1.1 franka\_interface

#### 1.1.1 ArmInterface

- Interface class that can monitor and control the robot
- Provides all required information about robot state and end-effector state
- Joint positions, velocities, and effort can be directly controlled and monitored using available methods
- Smooth interpolation of joint positions possible
- End-effector and Stiffness frames can be directly set (uses FrankaFramesInterface from franka\_ros\_interface/franka\_tools)

**class** franka\_interface.ArmInterface(synchronous\_pub=False)

Bases: `object`

Interface Class for an arm of Franka Panda robot Constructor.

**Parameters** `synchronous_pub` (`bool`) - designates the JointCommand Publisher as Synchronous if True and Asynchronous if False.

Synchronous Publishing means that all joint\_commands publishing to the robot's joints will block until the message has been serialized into a buffer and that buffer has been written to the transport of every current Subscriber. This yields predicable and consistent timing of messages being delivered from this Publisher. However, when using this mode, it is possible for a blocking Subscriber to prevent the joint\_command functions from exiting. Unless you need exact JointCommand timing, default to Asynchronous Publishing (False).

**class** RobotMode

Bases: `enum.IntEnum`

Enum class for specifying and retrieving the current robot mode.

**endpoint\_effort()**

Return Cartesian endpoint wrench {force, torque}.

**Return type** `dict({str:L{Limb.Point},str:L{Limb.Point}})`

**Returns**

force and torque at endpoint as named tuples in a dict

- 'force': Cartesian force on x,y,z axes in np.ndarray format
- 'torque': Torque around x,y,z axes in np.ndarray format

**endpoint\_pose()**

Return Cartesian endpoint pose {position, orientation}.

**Return type** `dict({str:L{Limb.Point},str:L{Limb.Quaternion}})`

**Returns**

position and orientation as named tuples in a dict

- 'position': np.array of x, y, z
- 'orientation': quaternion x,y,z,w in quaternion format

**endpoint\_velocity()**

Return Cartesian endpoint twist {linear, angular}.

**Return type** `dict({str:L{Limb.Point},str:L{Limb.Point}})`

**Returns**

linear and angular velocities as named tuples in a dict

- 'linear': np.array of x, y, z
- 'angular': np.array of x, y, z (angular velocity along the axes)

**error\_in\_current\_state()**

Return True if the specified limb has experienced an error.

**Return type** `bool`

**Returns** True if the arm has error, False otherwise.

**get\_joint\_limits()**

Return the joint limits (defined in the parameter server)

**Return type** `franka_core_msgs.msg.JointLimits`

**Returns** JointLimits

**get\_movegroup\_interface()**

**Returns** the movegroup interface instance associated with the robot.

**Return type** `franka_moveit.PandaMoveGroupInterface`

**get\_robot\_params()**

**Returns** Useful parameters from the ROS parameter server.

**Return type** `franka_interface.RobotParams`

**get\_robot\_status()**

Return dict with all robot status information.

**Return type** `dict`

**Returns** ['robot\_mode' (RobotMode object), 'robot\_status' (bool), 'errors' (dict() of errors and their truth value), 'error\_in\_curr\_status' (bool)]

**gravity\_comp()**

Return gravity compensation torques.

**Return type** np.ndarray

**Returns** 7D joint torques compensating for gravity.

**has\_collided()**

Returns true if either joint collision or cartesian collision is detected. Collision thresholds can be set using instance of [franka\\_tools.CollisionBehaviourInterface](#).

**in\_safe\_state()**

Return True if the specified limb is in safe state (no collision, reflex, errors etc.).

**Return type** bool

**Returns** True if the arm is in safe state, False otherwise.

**joint\_angle(joint)**

Return the requested joint angle.

**Parameters** **joint** (str) - name of a joint

**Return type** float

**Returns** angle in radians of individual joint

**joint\_angles()**

Return all joint angles.

**Return type** dict({str:float})

**Returns** unordered dict of joint name Keys to angle (rad) Values

**joint\_effort(joint)**

Return the requested joint effort.

**Parameters** **joint** (str) - name of a joint

**Return type** float

**Returns** effort in Nm of individual joint

**joint\_efforts()**

Return all joint efforts.

**Return type** dict({str:float})

**Returns** unordered dict of joint name Keys to effort (Nm) Values

**joint\_inertia\_matrix()**

Return joint inertia matrix (7,7)

**Return type** np.ndarray [7x7]

**joint\_names()**

Return the names of the joints for the specified limb.

**Return type** [str]

**Returns** ordered list of joint names from proximal to distal (i.e. shoulder to wrist).

**joint\_ordered\_angles()**

Return all joint angles.

**Return type** [double]

**Returns** joint angles (rad) orded by joint\_names from proximal to distal (i.e. shoulder to wrist).

**joint\_velocities()**

Return all joint velocities.

**Return type** dict({str:float})

**Returns** unordered dict of joint name Keys to velocity (rad/s) Values

**joint\_velocity(joint)**

Return the requested joint velocity.

**Parameters** joint (str) - name of a joint

**Return type** float

**Returns** velocity in radians/s of individual joint

**move\_to\_joint\_positions**(positions, timeout=10.0, threshold=0.00085,  
test=None, use\_moveit=True)

(Blocking) Commands the limb to the provided positions.

Waits until the reported joint state matches that specified.

This function uses a low-pass filter to smooth the movement.

**Parameters**

- **positions** (dict({str:float})) - joint\_name:angle command
- **timeout** (float) - seconds to wait for move to finish [15]
- **threshold** (float) - position threshold in radians across each joint when move is considered successful [0.008726646]
- **test** - optional function returning True if motion must be aborted
- **use\_moveit** (bool) - if set to True, and movegroup interface is available, move to the joint positions using moveit planner.

**move\_to\_neutral**(timeout=15.0, speed=0.15)

Command the Limb joints to a predefined set of “neutral” joint angles. From rosparam /franka\_control/neutral\_pose.

**Parameters**

- **timeout** (float) - seconds to wait for move to finish [15]
- **speed** (float) - ratio of maximum joint speed for execution default=0.15; range= [0.0-1.0]

**reset\_EE\_frame()**

Reset EE frame to default. (defined by FrankaFramesInterface.DEFAULT\_TRANSFORMATIONS.EE\_FRAME global variable defined above)

**Return type** [bool, str]

**Returns** [success status of service request, error msg if any]

**set\_EE\_frame(frame)**

Set new EE frame based on the transformation given by ‘frame’, which is the transformation matrix defining the new desired EE frame with respect to the flange frame. Motion controllers are stopped for switching



**Parameters** `frame` ([`float` (16,)] / `np.ndarray` (4x4)) - transformation matrix of new EE frame wrt flange frame (column major)

**Return type** [`bool`, `str`]

**Returns** [success status of service request, error msg if any]

**set\_ee\_frame\_to\_link**(`frame_name`, `timeout=5.0`)

Set new EE frame to the same frame as the link frame given by 'frame\_name' Motion controllers are stopped for switching

**Parameters** `frame_name` (`str`) - desired tf frame name in the tf tree

**Return type** [`bool`, `str`]

**Returns** [success status of service request, error msg if any]

**set\_collision\_threshold**(`cartesian_forces=None`, `joint_torques=None`)

Set Force Torque thresholds for deciding robot has collided.

**Returns** True if service call successful, False otherwise

**Return type** `bool`

**Parameters**

- **cartesian\_forces** ([`float`] size 6) - Cartesian force threshold for collision detection [x,y,z,R,P,Y] (robot motion stops if violated)
- **joint\_torques** ([`float`] size 7) - Joint torque threshold for collision (robot motion stops if violated)

**set\_command\_timeout**(`timeout`)

Set the timeout in seconds for the joint controller

**Parameters** `timeout` (`float`) - timeout in seconds

**set\_joint\_position\_speed**(`speed=0.3`)

Set ratio of max joint speed to use during joint position moves (only for `move_to_joint_positions`).

Set the proportion of maximum controllable velocity to use during joint position control execution. The default ratio is 0.3, and can be set anywhere from [0.0-1.0] (clipped). Once set, a speed ratio will persist until a new execution speed is set.

**Parameters** `speed` (`float`) - ratio of maximum joint speed for execution default= 0.3; range= [0.0-1.0]

**set\_joint\_positions**(`positions`)

Commands the joints of this limb to the specified positions.

**Parameters** `positions` ([`float`]) - ordered joint angles (from joint1 to joint7) to be commanded

**set\_joint\_positions\_velocities**(`positions`, `velocities`)

Commands the joints of this limb using specified positions and velocities using impedance control. Command at time  $t$  is computed as:

$$u_t = coriolis\_factor * coriolis\_t + K\_p * (positions - curr\_positions) + K\_d * (velocities - curr\_velocities)$$

**Parameters**

- **positions** ([`float`]) - desired joint positions as an ordered list corresponding to joints given by `self.joint_names()`

- **velocities** ([float]) - desired joint velocities as an ordered list corresponding to joints given by self.joint\_names()

**set\_joint\_torques**(torques)

Commands the joints of this limb to the specified torques.

**Parameters** **torques** (dict({str:float})) - joint\_name:torque command

**set\_joint\_velocities**(velocities)

Commands the joints of this limb to the specified velocities.

**Parameters** **velocities** (dict({str:float})) - joint\_name:velocity command

**tip\_states**()

Return Cartesian endpoint state for a given tip name

**Return type** TipState object

**Returns** pose, velocity, effort, effort\_in\_K\_frame

**what\_errors**()

Return list of error messages if there is error in robot state

**Return type** [str]

**Returns** list of names of current errors in robot state

**zero\_jacobian**()

Return end-effector jacobian (6,7)

**Return type** np.ndarray [6x7]

### 1.1.2 GripperInterface

- Interface class to monitor and control gripper
- Gripper open, close methods
- Grasp, move joints methods

```
class franka_interface.GripperInterface(gripper_joint_names=['panda_finger_joint1',
                                                             'panda_finger_joint2'],
                                         ns='franka_ros_interface',      cali-
                                         brate=False)
```

Bases: `object`

Interface class for the gripper on the Franka Panda robot.

#### Parameters

- **gripper\_joint\_names** ([str]) - Names of the finger joints
- **ns** (str) - base namespace of interface ('franka\_ros\_interface'/'panda\_simulator')
- **calibrate** (bool) - Attempts to calibrate the gripper when initializing class (defaults True)

**close**()

close gripper to till collision is detected. Note: This is not exactly doing what it should. The behaviour is faked by catching the error thrown when trying to grasp a very small object with a very small force. Since the gripper will actually hit the

object before it reaches the commanded width, we catch the feedback and send the gripper stop command to stop it where it is.

**Returns** True if command was successful, False otherwise.

**Return type** bool

**grasp**(width, force, speed=None, epsilon\_inner=0.005, epsilon\_outer=0.005, wait\_for\_result=True, cb=None)  
Grasps an object.

An object is considered grasped if the distance  $d$  between the gripper fingers satisfies  $(\text{ext}\{\text{width}\} - \text{ext}\{\text{epsilon\_inner}\}) < d < (\text{ext}\{\text{width}\} + \text{ext}\{\text{epsilon\_outer}\})$ .

**Parameters**

- **width** (float) - Size of the object to grasp. [m]
- **speed** (float) - Closing speed. [m/s]
- **force** (float) - Grasping force. [N]
- **epsilon\_inner** (float) - Maximum tolerated deviation when the actual grasped width is smaller than the commanded grasp width.
- **epsilon\_outer** (float) - Maximum tolerated deviation when the actual grasped width is wider than the commanded grasp width.
- **cb** - Optional callback function to use when the service call is done

**Returns** True if an object has been grasped, false otherwise.

**Return type** bool

**home\_joints**(wait\_for\_result=False)  
Performs homing of the gripper.

After changing the gripper fingers, a homing needs to be done. This is needed to estimate the maximum grasping width.

**Parameters** **wait\_for\_result** (bool) - if True, this method will block till response is recieved from server

**Returns** success

**Return type** bool

**joint\_effort**(joint)  
Return the requested joint effort.

**Parameters** **joint** (str) - name of a joint

**Return type** float

**Returns** effort in Nm of individual joint

**joint\_efforts**()  
Return all joint efforts.

**Return type** dict({str:float})

**Returns** unordered dict of joint name Keys to effort (Nm) Values

**joint\_names**()  
Return the names of the joints for the specified limb.

**Return type** [str]

**Returns** ordered list of joint names.

**joint\_ordered\_efforts()**

Return all joint efforts.

**Return type** [double]

**Returns** joint efforts ordered by joint\_names.

**joint\_ordered\_positions()**

Return all joint positions.

**Return type** [double]

**Returns** joint positions ordered by joint\_names.

**joint\_ordered\_velocities()**

Return all joint velocities.

**Return type** [double]

**Returns** joint velocities ordered by joint\_names.

**joint\_position(joint)**

Return the requested joint position.

**Parameters** **joint** (**str**) - name of a joint

**Return type** **float**

**Returns** position individual joint

**joint\_positions()**

Return all joint positions.

**Return type** **dict**({**str**:**float**})

**Returns** unordered dict of joint name Keys to pos

**joint\_velocities()**

Return all joint velocities.

**Return type** **dict**({**str**:**float**})

**Returns** unordered dict of joint name Keys to velocity (rad/s) Values

**joint\_velocity(joint)**

Return the requested joint velocity.

**Parameters** **joint** (**str**) - name of a joint

**Return type** **float**

**Returns** velocity in radians/s of individual joint

**move\_joints**(width, speed=None, wait\_for\_result=True)

Moves the gripper fingers to a specified width.

**Parameters**

- **width** (**float**) - Intended opening width. [m]
- **speed** (**float**) - Closing speed. [m/s]
- **wait\_for\_result** (**bool**) - if True, this method will block till response is recieved from server

**Returns** True if command was successful, False otherwise.

**Return type** `bool`

**open()**

Open gripper to max possible width.

**Returns** True if command was successful, False otherwise.

**Return type** `bool`

**set\_velocity(value)**

Set default value for gripper joint motions. Used for move and grasp commands.

**Parameters** **value** (`float`) - speed value [m/s]

**stop\_action()**

Stops a currently running gripper move or grasp.

**Returns** True if command was successful, False otherwise.

**Return type** `bool`

### 1.1.1.3 RobotEnable

- Interface class to reset robot when in recoverable error (use `enable_robot.py` script in scripts/)

**class** `franka_interface.RobotEnable`(`robot_params=None`)

Bases: `object`

Class `RobotEnable` - simple control/status wrapper around robot state

`enable()` - enable all joints `disable()` - disable all joints `reset()` - reset all joints, reset all jrcp faults, disable the robot `stop()` - stop the robot, similar to hitting the e-stop button

**Parameters** **robot\_params** (`RobotParams`) - A `RobotParams` instance (optional)

**disable()**

Disable all joints

**enable()**

Enable all joints

**state()**

Returns the last known robot state.

**Return type** `str`

**Returns** "Enabled"/"Disabled"

### 1.1.1.4 RobotParams

- Collects and stores all useful information about the robot from the ROS parameter server

**class** `franka_interface.RobotParams`

Bases: `object`

Interface class for essential ROS parameters on Intera robot.

**get\_joint\_names()**

Return the names of the joints for the specified limb from ROS parameter.

**Return type** `list[str]`

**Returns** ordered list of joint names from proximal to distal (i.e. shoulder to wrist). joint names for limb

**get\_robot\_name()**

Return the name of class of robot from ROS parameter.

**Return type** `str`

**Returns** name of the robot

## 1.2 franka\_moveit

### 1.2.1 PandaMoveGroupInterface

- Provides interface to control and plan motions using MoveIt in ROS.
- Simple methods to plan and execute joint trajectories and cartesian path.
- Provides easy reset and environment definition functionalities (See ExtendedPlanningSceneInterface below).

**class** `franka_moveit.PandaMoveGroupInterface`

**arm\_group**

**Returns** The MoveGroupCommander instance of this object. This is an interface to one group of joints. In this case the group is the joints in the Panda arm. This interface can be used to plan and execute motions on the Panda.

**Return type**

`moveit_commander.MoveGroupCommander`

---

**Note:** available\_methods: [http://docs.ros.org/jade/api/moveit\\_commander/html/classmoveit\\_\\_commander\\_1\\_1move\\_group\\_1\\_1MoveGroupCommander.html](http://docs.ros.org/jade/api/moveit_commander/html/classmoveit__commander_1_1move_group_1_1MoveGroupCommander.html)

---

**close\_gripper**(wait=False)

Using named states defined in urdf.

---

**Note:** If this named state is not found, your ros environment is probably not using the right panda\_moveit\_config package. Ensure that sourced package is from this repo -> [https://github.com/justagist/panda\\_moveit\\_config](https://github.com/justagist/panda_moveit_config)

---

**go\_to\_joint\_positions**(positions, wait=True, tolerance=0.005)

**Returns** status of joint motion plan execution

**Return type** `bool`

**Parameters**

- **positions** ([double]) - target joint positions (ordered)

- **wait** (bool) - if True, function will wait for trajectory execution to complete
- **tolerance** (double) - maximum error in final position for each joint to consider task a success

### **gripper\_group**

**Returns** The MoveGroupCommander instance of this object. This is an interface to one group of joints. In this case the group is the joints in the Panda arm. This interface can be used to plan and execute motions on the Panda.

#### **Return type**

moveit\_commander.MoveGroupCommander

---

**Note:** available\_methods: [http://docs.ros.org/jade/api/moveit\\_commander/html/classmoveit\\_commander\\_1\\_1move\\_group\\_1\\_1MoveGroupCommander.html](http://docs.ros.org/jade/api/moveit_commander/html/classmoveit_commander_1_1move_group_1_1MoveGroupCommander.html)

---

**move\_to\_neutral**(wait=True)

Send arm group to neutral pose defined using named state in urdf.

**open\_gripper**(wait=False)

Using named states defined in urdf.

---

**Note:** If this named state is not found, your ros environment is probably not using the right panda\_moveit\_config package. Ensure that sourced package is from this repo -> [https://github.com/justagist/panda\\_moveit\\_config](https://github.com/justagist/panda_moveit_config).

---

**plan\_joint\_path**(joint\_position)

:return plan for executing joint trajectory

### **robot\_state\_interface**

**Returns** The RobotCommander instance of this object

#### **Return type**

moveit\_commander.RobotCommander

---

**Note:** available\_methods: [http://docs.ros.org/jade/api/moveit\\_commander/html/classmoveit\\_commander\\_1\\_1robot\\_1\\_1RobotCommander.html](http://docs.ros.org/jade/api/moveit_commander/html/classmoveit_commander_1_1robot_1_1RobotCommander.html)

---

### **scene**

**Returns** The RobotCommander instance of this object. This is an interface to the world surrounding the robot

#### **Return type**

moveit\_commander.RobotCommander

---

**Note:** available\_methods: [http://docs.ros.org/indigo/api/moveit\\_ros\\_planning\\_interface/html/classmoveit\\_1\\_1planning\\_interface\\_1\\_1](http://docs.ros.org/indigo/api/moveit_ros_planning_interface/html/classmoveit_1_1planning_interface_1_1)

---

[1PlanningSceneInterface.html](#)

---

**set\_velocity\_scale**(value, group='arm')

Set the max velocity scale for executing planned motion. :param value: scale value (allowed (0,1] )

## 1.2.2 ExtendedPlanningSceneInterface

- Easily define scene for robot motion planning (MoveIt plans will avoid defined obstacles if possible).

**class** franka\_moveit.**ExtendedPlanningSceneInterface**

Bases: moveit\_commander.planning\_scene\_interface.PlanningSceneInterface

**add\_box**(name, pose, size, timeout=5)

Add object to scene and check if it is created.

**Parameters**

- **name** (**str**) - name of object
- **pose** (**geometry\_msgs.msg.PoseStamped**) - desired pose for the box
- **size** (**[float]** (len 3)) - size of the box
- **timeout** (**float**) - time in sec to wait while checking if box is created

**remove\_box**(box\_name, timeout=5)

Remove box from scene.

**Parameters**

- **box\_name** (**str**) - name of object
- **timeout** (**float**) - time in sec to wait while checking if box is created

## 1.3 franka\_tools

### 1.3.1 CollisionBehaviourInterface

- Define collision and contact thresholds for the robot safety and contact detection.

**class** franka\_tools.**CollisionBehaviourInterface**

Helper class to set collision and contact thresholds at cartesian and joint levels. (This class has no 'getter' functions to access the currently set collision behaviour values.)

**set\_collision\_threshold**(joint\_torques=None, cartesian\_forces=None)

**Returns** True if service call successful, False otherwise

**Return type** **bool**

**Parameters**

- **joint\_torques** (**[float]** size 7) - Joint torque threshold for collision (robot motion stops if violated)
- **cartesian\_forces** (**[float]** size 6) - Cartesian force threshold for collision detection [x,y,z,R,P,Y] (robot motion stops if violated)



**set\_contact\_threshold**(joint\_torques=None, cartesian\_forces=None)

**Returns** True if service call successful, False otherwise

**Return type** bool

**Parameters**

- **joint\_torques** ([float] size 7) - Joint torque threshold for identifying as contact
- **cartesian\_forces** ([float] size 6) - Cartesian force threshold for identifying as contact

**set\_force\_threshold\_for\_collision**(cartesian\_force\_values)

**Returns** True if service call successful, False otherwise

**Return type** bool

**Parameters** **cartesian\_force\_values** ([float] size 6) - Cartesian force threshold for collision detection [x,y,z,R,P,Y] (robot motion stops if violated)

**set\_force\_threshold\_for\_contact**(cartesian\_force\_values)

**Returns** True if service call successful, False otherwise

**Return type** bool

**Parameters** **cartesian\_force\_values** ([float] size 6) - Cartesian force threshold for contact detection [x,y,z,R,P,Y]

**set\_ft\_contact\_collision\_behaviour**(torque\_lower=None, torque\_upper=None, force\_lower=None, force\_upper=None)

**Returns** True if service call successful, False otherwise

**Return type** bool

**Parameters**

- **torque\_lower** ([float] size 7) - Joint torque threshold for contact detection
- **torque\_upper** ([float] size 7) - Joint torque threshold for collision (robot motion stops if violated)
- **force\_lower** ([float] size 6) - Cartesian force threshold for contact detection [x,y,z,R,P,Y]
- **force\_upper** ([float] size 6) - Cartesian force threshold for collision detection [x,y,z,R,P,Y] (robot motion stops if violated)

### 1.3.2 FrankaControllerManagerInterface

- List, start, stop, load available controllers for the robot
- Get the current controller status (commands, set points, controller gains, etc.)
- Update controller parameters through ControllerParamConfigClient (see below)

```
class franka_tools.FrankaControllerManagerInterface(ns='franka_ros_interface',
                                                    synchronous_pub=False,
                                                    sim=False)
```

Bases: `object`

### Parameters

- **synchronous\_pub** (`bool`) - designates the JointCommand Publisher as Synchronous if True and Asynchronous if False.

Synchronous Publishing means that all joint\_commands publishing to the robot's joints will block until the message has been serialized into a buffer and that buffer has been written to the transport of every current Subscriber. This yields predicable and consistent timing of messages being delivered from this Publisher. However, when using this mode, it is possible for a blocking Subscriber to prevent the joint\_command functions from exiting. Unless you need exact JointCommand timing, default to Asynchronous Publishing (False).

- **ns** (`str`) - base namespace of interface ('frank\_ros\_interface'/'panda\_simulator')
- **sim** (`bool`) - Flag specifying whether the robot is in simulation or not (can be obtained from `franka_interface.RobotParams` instance)

### `controller_dict()`

Get all controllers as dict

**Returns** name of the controller to be stopped

**Return type** dict {'controller\_name': ControllerState}

### `get_controller_config_client(controller_name)`

**Returns** The parameter configuration client object associated with the specified controller

**Return type** ControllerParamConfigClient obj (if None, returns False)

**Parameters controller\_name** (`str`) - name of controller whose config client is required

### `get_controller_state()`

Get the status of the current controller, including set points, computed command, controller gains etc. See the ControllerStateInfo class (above) parameters for more info.

### `get_current_controller_config_client()`

**Returns** The parameter configuration client object associated with the currently active controller

**Return type** ControllerParamConfigClient obj (if None, returns False)

**Parameters controller\_name** (`str`) - name of controller whose config client is required

### `is_loaded(controller_name)`

Check if the given controller is loaded.

**Parameters controller\_name** (`str`) - name of controller whose status is to be checked

**Returns** True if controller is loaded, False otherwise

**Return type** `bool`

**is\_running**(controller\_name)

Check if the given controller is running.

**Parameters** **controller\_name** (`str`) - name of controller whose status is to be checked

**Returns** True if controller is running, False otherwise

**Return type** `bool`

**list\_active\_controller\_names**(only\_motion\_controllers=False)

**Returns** List of names active controllers associated to a controller manager namespace.

**Return type** [`str`]

**Parameters** **only\_motion\_controller** (`bool`) - if True, only motion controllers are returned

**list\_active\_controllers**(only\_motion\_controllers=False)

**Returns** List of active controllers associated to a controller manager namespace. Contains both stopped/running controllers, as returned by the `list_controllers` service, plus uninitialized controllers with configurations loaded in the parameter server.

**Return type** [ControllerState obj]

**Parameters** **only\_motion\_controller** (`bool`) - if True, only motion controllers are returned

**list\_controller\_names**()

**Returns** List of names all controllers associated to a controller manager namespace.

**Return type** [`str`]

**Parameters** **only\_motion\_controller** (`bool`) - if True, only motion controllers are returned

**list\_controller\_types**()

**Returns** List of controller types associated to a controller manager namespace. Contains both stopped/running/loaded controllers, as returned by the `list_controller_types` service, plus uninitialized controllers with configurations loaded in the parameter server.

**Return type** [`str`]

**list\_controllers**()

**Returns** List of controllers associated to a controller manager namespace. Contains both stopped/running controllers, as returned by the `list_controllers` service, plus uninitialized controllers with configurations loaded in the parameter server.

**Return type** [ControllerState obj]

**list\_loaded\_controllers**()

**Returns** List of controller types associated to a controller manager namespace. Contains all loaded controllers, as returned by the `list_controller_types` service, plus uninitialized controllers with configurations loaded in the parameter server.

**Return type** `[str]`

**`list_motion_controllers()`**

**Returns** List of motion controllers associated to a controller manager namespace. Contains both stopped/running controllers, as returned by the `list_controllers` service, plus uninitialized controllers with configurations loaded in the parameter server.

**Return type** `[ControllerState obj]`

**`load_controller(name)`**

Loads the specified controller

**Parameters** `name (str)` - name of the controller to be loaded

**`set_motion_controller(controller_name)`**

Set the specified controller as the (only) motion controller

**Returns** name of currently active controller (can be used to switch back to this later)

**Return type** `str`

**Parameters** `controller_name (str)` - name of controller to start

**`start_controller(name)`**

Starts the specified controller

**Parameters** `name (str)` - name of the controller to be started

**`stop_controller(name)`**

Stops the specified controller

**Parameters** `name (str)` - name of the controller to be stopped

**`unload_controller(name)`**

Unloads the specified controller

**Parameters** `name (str)` - name of the controller to be unloaded

### 1.3.3 ControllerParamConfigClient

- Get and set the controller parameters (gains) for the active controller

**`class franka_tools.ControllerParamConfigClient(controller_name)`**

Interface class for updating dynamically configurable parameters of a controller.

**Parameters** `controller_name (str)` - The name of the controller.

**`get_config(timeout=5)`**

**Returns** the currently set values for all parameters from the server

**Return type** `dict {str : float}`

**Parameters** `timeout (float)` - time to wait before giving up on service request

**get\_controller\_gains**(timeout=5)

**Returns** the currently set values for controller gains from the server

**Return type** ( [float], [float] )

**Parameters** **timeout** (float) - time to wait before giving up on service request

**get\_joint\_motion\_smoothing\_parameter**(timeout=5)

**Returns** the currently set value for the joint position smoothing parameter from the server.

**Return type** float

**Parameters** **timeout** (float) - time to wait before giving up on service request

**get\_parameter\_descriptions**(timeout=5)

**Returns** the description of each parameter as defined in the cfg file from the server.

**Return type** dict {str : str}

**Parameters** **timeout** (float) - time to wait before giving up on service request

**is\_running**

**Returns** True if client is running / server is unavailable; False otherwise

**Return type** bool

**set\_controller\_gains**(k\_gains, d\_gains=None)

Update the stiffness and damping parameters of the joints for the current controller.

**Parameters**

- **k\_gains** ([float]) - joint stiffness parameters (should be within limits specified in franka documentation; same is also set in franka\_ros\_controllers/cfg/joint\_controller\_params.cfg)
- **d\_gains** ([float]) - joint damping parameters (should be within limits specified in franka documentation; same is also set in franka\_ros\_controllers/cfg/joint\_controller\_params.cfg)

**set\_joint\_motion\_smoothing\_parameter**(value)

**Update the joint motion smoothing parameter (only valid for position\_joint\_position\_controller).**

**Parameters** **value** ([float]) - smoothing factor (should be within limit set in franka\_ros\_controllers/cfg/joint\_controller\_params.cfg)

**start**(timeout=5)

Start the dynamic\_reconfigure client

**Parameters** **timeout** (float) - time to wait before giving up on service request

**update\_config**(\*\*kwargs)

Update the config in the server using the provided keyword arguments.

**Parameters kwargs** – These are keyword arguments matching the parameter names in config file: franka\_ros\_controllers/cfg/joint\_controller\_params.cfg

### 1.3.4 FrankaFramesInterface

- Get and Set end-effector frame and stiffness frame of the robot easily
- Set the frames to known frames (such as links on the robot) directly

#### **class** franka\_tools.FrankaFramesInterface

Helper class to retrieve and set EE frames

Has to be updated externally each time franka states is updated. This is done by default within the PandaArm class (panda\_robot package: [https://github.com/justagist/panda\\_robot](https://github.com/justagist/panda_robot)).

Note that all controllers have to be unloaded before switching frames. This has to be done externally (also automatically handled in PandaArm class).

#### **frames\_are\_same**(frame1, frame2)

@return True if two transformation matrices are equal @rtype: bool @param frame1: 4x4 transformation matrix representing frame1 @type frame1: np.ndarray (shape 4x4), or list (flattened column major 4x4) @param frame2: 4x4 transformation matrix representing frame2 @type frame2: np.ndarray (shape 4x4), or list (flattened column major 4x4)

#### **get\_EE\_frame**(as\_mat=False)

Get current EE frame transformation matrix in flange frame

@type as\_mat: bool @param as\_mat: if True, return np array, else as list @rtype: [float (16,)] / np.ndarray (4x4) @return: transformation matrix of EE frame wrt flange frame (column major)

#### **get\_K\_frame**(as\_mat=False)

Get current K frame transformation matrix in EE frame

@type as\_mat: bool @param as\_mat: if True, return np array, else as list @rtype: [float (16,)] / np.ndarray (4x4) @return: transformation matrix of K frame wrt EE frame

#### **reset\_EE\_frame**()

Reset EE frame to default. (defined by DEFAULT\_TRANSFORMATIONS.EE\_FRAME global variable defined above)

@rtype: bool @return: success status of service request

#### **reset\_K\_frame**()

Reset K frame to default. (defined by **DEFAULT\_K\_FRAME** global variable defined above)

@rtype: bool @return: success status of service request

#### **set\_EE\_frame**(frame)

Set new EE frame based on the transformation given by 'frame', which is the transformation matrix defining the new desired EE frame with respect to the flange frame.

@type frame: [float (16,)] / np.ndarray (4x4) @param frame: transformation matrix of new EE frame wrt flange frame (column major) @rtype: bool @return: success status of service request

**set\_EE\_frame\_to\_link**(frame\_name, timeout=5.0)

Set new EE frame to the same frame as the link frame given by 'frame\_name' Motion controllers are stopped for switching

@type frame\_name: str @param frame\_name: desired tf frame name in the tf tree  
@rtype: [bool, str] @return: [success status of service request, error msg if any]

**set\_K\_frame**(frame)

Set new K frame based on the transformation given by 'frame', which is the transformation matrix defining the new desired K frame with respect to the EE frame.

@type frame: [float (16,)] / np.ndarray (4x4) @param frame: transformation matrix of new K frame wrt EE frame @rtype: bool @return: success status of service request

**set\_K\_frame\_to\_link**(frame\_name, timeout=5.0)

Set new K frame to the same frame as the link frame given by 'frame\_name' Motion controllers are stopped for switching

@type frame\_name: str @param frame\_name: desired tf frame name in the tf tree  
@rtype: [bool, str] @return: [success status of service request, error msg if any]

### 1.3.5 JointTrajectoryActionClient

- Command robot to given joint position(s) smoothly. (Uses the FollowJointTrajectory service from ROS control\_msgs package)
- Smoothly move to a desired (valid) pose without having to interpolate for smoothness (trajectory interpolation done internally)

```
class franka_tools.JointTrajectoryActionClient(joint_names,  
                                              ns='franka_ros_interface', con-  
                                              troller_name='position_joint_trajectory_controller')
```

Bases: `object`





## Indices and tables

- `genindex`
- `modindex`
- `search`



## Python Module Index

### f

franka\_interface, [1](#)  
franka\_moveit, [10](#)  
franka\_tools, [12](#)



## A

add\_box()  
     (franka\_moveit.ExtendedPlanningSceneInterface  
     method), 12  
 arm\_group (franka\_moveit.PandaMoveGroupInterface  
     attribute), 10  
 ArmInterface (class in franka\_interface), 1  
 ArmInterface.RobotMode (class in franka\_interface), 1

## C

close() (franka\_interface.GripperInterface method), 6  
 close\_gripper()  
     (franka\_moveit.PandaMoveGroupInterface  
     method), 10  
 CollisionBehaviourInterface (class in franka\_tools),  
     12  
 controller\_dict()  
     (franka\_tools.FrankaControllerManagerInterface  
     method), 14  
 ControllerParamConfigClient (class in franka\_tools),  
     16

## D

disable() (franka\_interface.RobotEnable method), 9

## E

enable() (franka\_interface.RobotEnable method), 9  
 endpoint\_effort() (franka\_interface.ArmInterface  
     method), 1  
 endpoint\_pose() (franka\_interface.ArmInterface  
     method), 2  
 endpoint\_velocity() (franka\_interface.ArmInterface  
     method), 2  
 error\_in\_current\_state()  
     (franka\_interface.ArmInterface method), 2  
 ExtendedPlanningSceneInterface (class in  
     franka\_moveit), 12

## F

frames\_are\_same()  
     (franka\_tools.FrankaFramesInterface  
     method), 18  
 franka\_interface (module), 1  
 franka\_moveit (module), 10

franka\_tools (module), 12  
 FrankaControllerManagerInterface (class in  
     franka\_tools), 13  
 FrankaFramesInterface (class in franka\_tools), 18

## G

get\_config()  
     (franka\_tools.ControllerParamConfigClient  
     method), 16  
 get\_controller\_config\_client()  
     (franka\_tools.FrankaControllerManagerInterface  
     method), 14  
 get\_controller\_gains()  
     (franka\_tools.ControllerParamConfigClient  
     method), 16  
 get\_controller\_state()  
     (franka\_tools.FrankaControllerManagerInterface  
     method), 14  
 get\_current\_controller\_config\_client()  
     (franka\_tools.FrankaControllerManagerInterface  
     method), 14  
 get\_EE\_frame() (franka\_tools.FrankaFramesInterface  
     method), 18  
 get\_joint\_limits() (franka\_interface.ArmInterface  
     method), 2  
 get\_joint\_motion\_smoothing\_parameter()  
     (franka\_tools.ControllerParamConfigClient  
     method), 17  
 get\_joint\_names() (franka\_interface.RobotParams  
     method), 9  
 get\_K\_frame() (franka\_tools.FrankaFramesInterface  
     method), 18  
 get\_movegroup\_interface()  
     (franka\_interface.ArmInterface method), 2  
 get\_parameter\_descriptions()  
     (franka\_tools.ControllerParamConfigClient  
     method), 17  
 get\_robot\_name() (franka\_interface.RobotParams  
     method), 10  
 get\_robot\_params() (franka\_interface.ArmInterface  
     method), 2  
 get\_robot\_status() (franka\_interface.ArmInterface  
     method), 2  
 go\_to\_joint\_positions()  
     (franka\_moveit.PandaMoveGroupInterface  
     method), 10

grasp() (franka\_interface.GripperInterface method), 7  
 gravity\_comp() (franka\_interface.ArmInterface method), 2  
 gripper\_group (franka\_moveit.PandaMoveGroupInterface attribute), 11  
 GripperInterface (class in franka\_interface), 6

## H

has\_collided() (franka\_interface.ArmInterface method), 3  
 home\_joints() (franka\_interface.GripperInterface method), 7

## I

in\_safe\_state() (franka\_interface.ArmInterface method), 3  
 is\_loaded() (franka\_tools.FrankaControllerManagerInterface method), 14  
 is\_running (franka\_tools.ControllerParamConfigClient attribute), 17  
 is\_running() (franka\_tools.FrankaControllerManagerInterface method), 15

## J

joint\_angle() (franka\_interface.ArmInterface method), 3  
 joint\_angles() (franka\_interface.ArmInterface method), 3  
 joint\_effort() (franka\_interface.ArmInterface method), 3  
 joint\_effort() (franka\_interface.GripperInterface method), 7  
 joint\_efforts() (franka\_interface.ArmInterface method), 3  
 joint\_efforts() (franka\_interface.GripperInterface method), 7  
 joint\_inertia\_matrix() (franka\_interface.ArmInterface method), 3  
 joint\_names() (franka\_interface.ArmInterface method), 3  
 joint\_names() (franka\_interface.GripperInterface method), 7  
 joint\_ordered\_angles() (franka\_interface.ArmInterface method), 3  
 joint\_ordered\_efforts() (franka\_interface.GripperInterface method), 8  
 joint\_ordered\_positions() (franka\_interface.GripperInterface method), 8  
 joint\_ordered\_velocities() (franka\_interface.GripperInterface method), 8  
 joint\_position() (franka\_interface.GripperInterface method), 8  
 joint\_positions() (franka\_interface.GripperInterface method), 8  
 joint\_velocities() (franka\_interface.ArmInterface method), 4  
 joint\_velocities() (franka\_interface.GripperInterface method), 8  
 joint\_velocity() (franka\_interface.ArmInterface method), 4  
 joint\_velocity() (franka\_interface.GripperInterface method), 8

JointTrajectoryActionClient (class in franka\_tools), 19

## L

list\_active\_controller\_names() (franka\_tools.FrankaControllerManagerInterface method), 15  
 list\_active\_controllers() (franka\_tools.FrankaControllerManagerInterface method), 15  
 list\_controller\_names() (franka\_tools.FrankaControllerManagerInterface method), 15  
 list\_controller\_types() (franka\_tools.FrankaControllerManagerInterface method), 15  
 list\_controllers() (franka\_tools.FrankaControllerManagerInterface method), 15  
 list\_loaded\_controllers() (franka\_tools.FrankaControllerManagerInterface method), 15  
 list\_motion\_controllers() (franka\_tools.FrankaControllerManagerInterface method), 16  
 load\_controller() (franka\_tools.FrankaControllerManagerInterface method), 16

## M

move\_joints() (franka\_interface.GripperInterface method), 8  
 move\_to\_joint\_positions() (franka\_interface.ArmInterface method), 4  
 move\_to\_neutral() (franka\_interface.ArmInterface method), 4  
 move\_to\_neutral() (franka\_moveit.PandaMoveGroupInterface method), 11

## O

open() (franka\_interface.GripperInterface method), 9  
 open\_gripper() (franka\_moveit.PandaMoveGroupInterface method), 11

## P

PandaMoveGroupInterface (class in franka\_moveit), 10  
 plan\_joint\_path() (franka\_moveit.PandaMoveGroupInterface method), 11

## R

remove\_box() (franka\_moveit.ExtendedPlanningSceneInterface method), 12  
 reset\_EE\_frame() (franka\_interface.ArmInterface method), 4  
 reset\_EE\_frame() (franka\_tools.FrankaFramesInterface method), 18  
 reset\_K\_frame() (franka\_tools.FrankaFramesInterface method), 18

robot\_state\_interface  
     (franka\_moveit.PandaMoveGroupInterface  
     attribute), 11  
 RobotEnable (class in franka\_interface), 9  
 RobotParams (class in franka\_interface), 9

## S

scene (franka\_moveit.PandaMoveGroupInterface  
     attribute), 11  
 set\_collision\_threshold()  
     (franka\_interface.ArmInterface method), 5  
 set\_collision\_threshold()  
     (franka\_tools.CollisionBehaviourInterface  
     method), 12  
 set\_command\_timeout()  
     (franka\_interface.ArmInterface method), 5  
 set\_contact\_threshold()  
     (franka\_tools.CollisionBehaviourInterface  
     method), 13  
 set\_controller\_gains()  
     (franka\_tools.ControllerParamConfigClient  
     method), 17  
 set\_EE\_frame() (franka\_interface.ArmInterface  
     method), 4  
 set\_EE\_frame() (franka\_tools.FrankaFramesInterface  
     method), 18  
 set\_EE\_frame\_to\_link()  
     (franka\_interface.ArmInterface method), 5  
 set\_EE\_frame\_to\_link()  
     (franka\_tools.FrankaFramesInterface  
     method), 18  
 set\_force\_threshold\_for\_collision()  
     (franka\_tools.CollisionBehaviourInterface  
     method), 13  
 set\_force\_threshold\_for\_contact()  
     (franka\_tools.CollisionBehaviourInterface  
     method), 13  
 set\_ft\_contact\_collision\_behaviour()  
     (franka\_tools.CollisionBehaviourInterface  
     method), 13  
 set\_joint\_motion\_smoothing\_parameter()  
     (franka\_tools.ControllerParamConfigClient  
     method), 17  
 set\_joint\_position\_speed()  
     (franka\_interface.ArmInterface method), 5  
 set\_joint\_positions()  
     (franka\_interface.ArmInterface method), 5  
 set\_joint\_positions\_velocities()  
     (franka\_interface.ArmInterface method), 5  
 set\_joint\_torques() (franka\_interface.ArmInterface  
     method), 6  
 set\_joint\_velocities()  
     (franka\_interface.ArmInterface method), 6  
 set\_K\_frame() (franka\_tools.FrankaFramesInterface  
     method), 19  
 set\_K\_frame\_to\_link()  
     (franka\_tools.FrankaFramesInterface  
     method), 19  
 set\_motion\_controller()  
     (franka\_tools.FrankaControllerManagerInterface  
     method), 16  
 set\_velocity() (franka\_interface.GripperInterface  
     method), 9  
 set\_velocity\_scale()  
     (franka\_moveit.PandaMoveGroupInterface  
     method), 12

start() (franka\_tools.ControllerParamConfigClient  
     method), 17  
 start\_controller()  
     (franka\_tools.FrankaControllerManagerInterface  
     method), 16  
 state() (franka\_interface.RobotEnable method), 9  
 stop\_action() (franka\_interface.GripperInterface  
     method), 9  
 stop\_controller()  
     (franka\_tools.FrankaControllerManagerInterface  
     method), 16

## T

tip\_states() (franka\_interface.ArmInterface method),  
     6

## U

unload\_controller()  
     (franka\_tools.FrankaControllerManagerInterface  
     method), 16  
 update\_config()  
     (franka\_tools.ControllerParamConfigClient  
     method), 17

## W

what\_errors() (franka\_interface.ArmInterface  
     method), 6

## Z

zero\_jacobian() (franka\_interface.ArmInterface  
     method), 6